FT-IR Spectrometer - Nicolet iS10 Department of Chemistry Pondicherry University, Puducherry - 605014

SAMPLE ANALYSIS REQUISITION FORM

Name of the Student/Scholar

Signature of the Faculty - In charge:

Course Registered	:		
Department / Institution	:		
Name of the Supervisor/Guide	:		
Signature of the Supervisor	:		Date:
Signature of the Head of the Department	:		
Number of Samples	:		
Sample Code	:		
Nature of Sample	:	Solid/ Liquid/Solution	
Solvent information	:		
Note:			
(i) Entry of moisture containing materials are strictly prohibited inside FT - IR room.			
(ii) Dry Solvents must be used for solution cell usage since KBr cells are sensitive to moisture.			